


ITEM NO.	QTY.	PART NUMBER	DESCRIPTION	MATERIAL	Number
1	2	15.2925	NIPPLE .25 HP X 02.25 316 SS	316 SS	ND
2	1	15.2400	TEE .25 HP OLDSTYL .75THK 316 SS	316 SS	ND
3	1	1.3700	GAUGE, 0-35000 MCDANIEL	316 SS	A-10554
4	1	40.10359	TUBE BENT .12 HP GAUGE OPTION LM20	316 SS	40.10359
5	2	15.2000	GLAND .25 HP	316 SS	ND
6	2	15.2100	COLLAR.25 HP	316 SS	ND
7	1	25.6019	ADPT .25 HPM X .12 HPF	NITRONIC 60	B-12933
8	1	15.1900	GLAND, .13 HP	316 SS	ND
9	1	15.1901	COLLAR, .13 HP	316 SS	ND
10	1	REF	CHAMBER LV1	SS	ND

REVISIONS						
REV.	ZONE	ECN NO.	DESCRIPTION	DRAWN BY	APPROVED BY	DATE
A	-	-	RELEASED FOR MANUFACTURE	JR	-	10/6/15



MICROFLUIDICS
 90 GLACIER DRIVE
 SUITE 1000
 WESTWOOD, MA 02090 USA
 SCALE: 1:1
 DATE: 10/6/15
 TOLERANCES: UNLESS OTHERWISE NOTED
 = ± 1/2" XX = ± 1/32" XXX = ± .005"

PRESSURE GAUGE OPTION LM20
 MATERIAL: SEE BILL OF MATERIALS
 SIZE: B
 DRAWING NUMBER: 90.10111
 REV: A